

SYSTEMZ

 **TOPRO**





Made in Norway

THE NEW GENERATION OF ROLLATORS!

Increased safety, flexibility and the best possible comfort

**TOPRO
TROJA 2G**

The model on this page is TOPRO TROJA 2G Premium.



1 TOPRO Ergo Grip

Unique ergonomic grip that makes it easier to stand from a sitting position. Provides a comfortable armrest and offers different holding positions when standing and resting.

2 TOPRO Quick Release

Easy to exchange or take off the wheels. Choose between:

3 Comfort Wheels Soft (PUR)

Comfortable and secure on un-surfaced roads, cobblestones etc.

4 Comfort Wheels Grip (TPE)

Maneuverable with good grip on smooth surfaces like flooring and tarmac.

5 TOPRO Bell

Attract attention when necessary. Easy to use without having to take your hands off the handles.

6 TOPRO Shopping Bag or Net Basket

Both hold up to 5 kg and are easy to mount and remove. The Shopping bag includes a shoulder strap for easy handling. It has two pockets inside for small personal belongings and a lid to protect your items.

7 TOPRO LED light and fall alarm

Feel safe when walking in dark or low light conditions. The integrated motion and light sensors turn the lights on automatically. If the rollator should fall the lights will flash and emit a warning tone.

8 Reflective Panels

Be visible from all sides in the dark. The 3M Scotchlite Reflective panels increase safety.

In addition are the original features kept:

- Light weight
- Memory function for easy height adjustment of the handles.
- Tilt function and edge guard for safe and easy manoeuvring.
- Easy to fold.
- Stays firmly on its wheels when folded and takes little space.
- Tested and approved according to EN ISO 11199
- Biocompatibility tested according to EN ISO 10993-5

Choose from three different models

PREMIUM PLUS • Shopping bag • Comfort Wheels Soft (PUR) • LED light with alarm
PREMIUM • Shopping bag • Comfort Wheels Soft (PUR)
BASIC • Net basket • Comfort Wheels Grip (TPE)

The rest of the features are standard on all models

**TOPRO
TROJA 2G**

The model on this page is TOPRO TROJA 2G Premium Plus.

The innovative new TOPRO TROJA 2G stands ahead of the competition with increased safety, flexibility and the best possible comfort.

Its predecessor, the TOPRO TROJA Classic, was selected as winner of the highly respected German consumer test 'Stiftung Warentest (edition 9/2005)' since then the TOPRO engineers have been working to develop the rollator based on research and feedback from users and therapists. The result is the second generation TOPRO TROJA with even more user friendly functions, a high focus on safety, flexibility and comfort.

Statements from professionals about the TOPRO TROJA 2G Ergo Grip:



Fabian Crüwell, physiotherapist and executive therapist
Specialist Hospital Neresheim, Germany

"It gives a number of grip options. The grip requires more tension in the upper body, which is good for the patient. Rising up is easier, because you can support yourself better with this grip."



Yasko Hupfauer, physiotherapist
University and Rehabilitation Clinic Ulm, Germany

"Patients with deformities in their hands can hold the rollator easier. The user walks WITHIN the rollator, not behind."



Nora Kern, NOI and Bobath instructor, Germany
"Sitting on the walker is more comfortable, since the arms can easily rest on the Ergo Grip."

TOPRO



Size	M (medium)		S (small)	
	KG	LBS	KG	LBS
Max. user weight	150	330	125	275
Weight with basket and wheels (TPE)	7,7	17	7,3	16
Transport weight (without wheels and basket)	6	13	5,6	12
	CM	INCHES	CM	INCHES
Recommended user height	150 - 200	4ft 11" - 6ft 6 3/4"	135 - 170	4ft 5" - 5ft 7"
Height of support points	78 - 100	30 3/4" - 39 3/8"	67 - 86	26 3/8" - 33 3/4"
Maximum length	65	25 1/2"	65	25 1/2"
Maximum width	61	24"	61	24"
Distance between support points	47	18 1/2"	47	18 1/2"
Seat height	62	24 1/2"	54	21 1/4"
Turning diameter	84	33"	84	33"
Diameter front/rear wheels	20	8"	20	8"
Width front/rear wheels (TPE wheels)	3,6	1 1/2"	3,6	1 1/2"
Distance driving brake	6,5	2 1/2"	6,5	2 1/2"
Length when folded	65	26"	65	26"
Width when folded	24	9 1/2"	24	9 1/2"
Height when folded	81,5	32 1/2"	74	29 3/4"

Model / Colour		Dark grey	Dark blue*	Sand*		Silver
PREMIUM PLUS	Art no. / Model no.	814600 / 330	335	340	814601 /	430
PREMIUM	Art no. / Model no.	814600 / 310	315	320	814601 /	400
BASIC	Art no. / Model no.	814600 / 300			814601 /	415

Imperial measurements are approximate. * Limited edition colours. The colours may vary from the colours of this print.

Accessories for a variety of uses

To be ordered separately.

Product TOPRO ART. No.



Adjustable back support
with padding 814607



Tray 814728



Anti slip mat
for tray 814059



Crutch holder 814025



Bottle holder 814043



Transport bag 814042



Comfort Wheels Soft (PUR) 814610
Set of 4



Comfort Wheels Grip (TPE) 814617
Set of 4



Shopping bag with
cover and shoulder strap 814622



Net basket 814625



Rear bag with zipper 814046



LED light with
integrated fall alarm 814609



Lamp 814608



Nameplate 814024



Ice spurs for shoes
Size M (shoe size 36-40) 814029
Size L (shoe size 41-46) 814035



Umbrella 814735
(not for sale in Germany)



Lock 814023



Guiding handle for
the visually impaired 814014



Basket for oxygen bottle 814009
(only compressed, not liquid)



One-handed brake 814026



Forearm support (pair)
for TROJA M 814829
for TROJA S and TROJA X 814845

SYSTMZ

 **TOPRO**



Made in Norway

THE ORIGINAL!

A supreme test winner on the market since 2002

The acknowledged German test institute, Stiftung Warentest, tested 14 rollators in edition 9/2005. This is what they said about the Norwegian manufactured test winner, TOPRO TROJA:

- Best operating comfort.
- Easy to fold.
- Good brake function.
- Light weight and easy to manoeuvre.
- Requires little space when folded.
- The test team's favorite.

**TOPRO
TROJA Classic**



TOPRO TROJA ROLLATORS are easy to fold, requiring little space when stored

- Pull the cord on the seat and the rollator folds together horizontally.
- Push the handlebars down into the sidebars.
- Unique memory function that sets the handles back to the right height.
- Stands securely on its wheels when folded.
- Takes little space and is very practical to transport.
- Excellent for users who travel by car, train, bus or airplane.

An elegant and practical light weight rollator

- Plastic coated aluminium tube and plastic couplings.
- A total of 12 ball bearings in the wheels and forks provide excellent driving comfort.
- Can be hosed and cleaned.
- All wear and tear parts can be replaced and brakes can be adjusted.
- Tested and approved according to EN ISO 11199

Excellent brakes, kerb lifter and edge guard

- Easy to negotiate thresholds and kerbs etc.
- Safe to use.

Three sizes!

Colour

Maximum weight of user

Product weight

Recommended user height

Height of handles

Maximum length

Maximum width

Width between handles

Seat height

Turning diameter

Diameter front/rear wheels

Width front/rear wheels

Grip distance driving brake

Height when folded up

Width when folded up

Length when folded up

Handles

Area of use

**TOPRO
TROJA Classic**

TOPRO ART. No.



M (medium)

Dark grey

KG	LBS
165	364
7,4	16,3
CM	INCHES
150-200	4ft 11"-6ft 4 3/8"
78-100	30 3/4"-39 3/8"
65	25 5/8"
61	24
47	18 1/2"
62	24 3/8"
84	33 1/8"
20	7 7/8"
3,6	1 3/8"
6,5	2 5/8"
80	31 1/2"
24	9 1/2"
66	26

Anatomic

Designed and approved for both indoor and outdoor use

814750

S (small)

Silver

KG	LBS
125	274
7,1	15,6
CM	INCHES
135-170	4ft 5 1/8"-5ft 7"
67-86	26 3/8"-33 7/8"
65	25 5/8"
61	24
47	18 1/2"
54	21 1/4"
84	33 1/8"
20	7 7/8"
3,6	1 3/8"
6,5	2 5/8"
74	29 1/8"
24	9 1/2"
66	26

Anatomic

814770

X (extra small)

Silver

KG	LBS
80	176
6,6	14,2
CM	INCHES
125-160	4ft 1 1/2"-5ft 3"
63,5-79	25"-31"
65	25 5/8"
51,5	20 1/4"
37	14 1/2"
52	20 1/2"
84	33 1/8"
20	7 7/8"
3,6	1 3/8"
6,5	2 5/8"
65	25 5/8"
25	9 3/8"
66	26

Anatomic

814740



Also available complete with forearm supports. See TOPRO Walkers brochure.

Accessories for a variety of uses

To be ordered separately.

Product **TOPRO ART. No.**



Adjustable back support **814726**
with padding, for rollators from 2009



Back Support for
TROJA Walker **814765**



Tray **X** **814728**



Anti slip mat
for tray **X** **814059**



Crutch holder **814025**



Bottle holder **814043**



Shopping bag with
cover and shoulder strap **X** **814140**



Rear bag
with zipper **X** **814046**



Transport bag **814042**



Flower stickers,
set with 4 stickers **814141**



Reflective Panels **814624**



Lamp **814608**



Nameplate **814024**



Ice spurs for shoes
Size M (shoe size 36-40) **814029**
Size L (shoe size 41-46) **814035**



Umbrella **814735**
(not for sale in Germany) **157207**



Lock **814023**



Guiding handle for
the visually impaired **814014**



Basket for oxygen bottle **814009**
(only compressed, not liquid) **X**



One-handed brake **814026**
157197



Forearm support (pair)
for TROJA M **814829**
for TROJA S and TROJA X **814845**

Singapore (HQ)

Systmz Pte Ltd
No 1, Bukit Batok Crescent
WCEGA Plaza #09-07
Singapore 658064
Tel: +65 6465 0387
Fax: +65 6465 0656
Email: info@systmz.sg

Malaysia

Systmz Sdn Bhd
No.11, Jalan PJU 3/46
Sunway Damansara Technology Park
47810 Petaling Jaya
Selangor Darul Ehsan
Malaysia
Customer Care: 1300 222 500
Email: info@systmz.com.my

Dubai Liason Office

Al Leddan General Trading
P.O. Box 488176,
Dubai - UAE, Barsha
Dusseldorf Business Tower - 402
Tel: +971 4 4570375
Mobile: +971 50 2323001
Email: gm@leddan-ics.com

SYSTEMZ

 **TOPRO**



STEADY!

Safe walkers for institutions and at home

**TOPRO
TAURUS**

Steady and safe walkers for institutions and at home

- Easy to read heightscales and easy to adjust the height.
- High enough for tall users (up to 2 m).
- Tested and approved according to ISO 11199-3 for use indoors on plain foundation.
- The TOPRO TAURUS E-models are tested and approved according to IEC 60601-1:2005 + corr 1:2006 + corr 2:2007 and IEC 60601-1-2:2007

Maintenance free material

- Aluminium with polyester surface and polyurethane.
- Easy to clean.
- Non-corrosive.
- All wear and tear parts can be exchanged.

Two height adjustment methods:

TOPRO TAURUS H Hydraulic height adjustment.

TOPRO TAURUS E Electrical height adjustment.

Two levels of equipment:

BASIC Standard equipment with adjustable arm pads and foot brakes on all wheels. Delivered partly assembled.

PREMIUM Premium has in addition:



		TOPRO ART. No.
TOPRO TAURUS H	BASIC	814785
TOPRO TAURUS H	PREMIUM	814787
TOPRO TAURUS E	BASIC	814790
TOPRO TAURUS E	PREMIUM	814789

TOPRO TAURUS H

- Gas springs simplify the height adjustment.
- Can be used in the shower.

Product weight 18 kg / 2 ¾ Stone

Max user weight 150 kg / 23 ½ Stone

Imperial measurements are approximate.

**TOPRO
TAURUS**





TOPRO TAURUS E

- Electric motor makes the height adjustment very simple.
- Can also be used as raising support.

Product weight 22 kg / 3 1/2 Stone

Max user weight 150 kg / 23 1/2 Stone

Imperial measurements are approximate.

Accessories for a variety of uses

To be ordered separately.

Product	TOPRO ART. No.			
			Foot plate	814645
Tray	814810		Height adjustment panel (Concens)	814646
			Wall mount for the Taurus battery charger (Concens)	814085
Brake handle left side	814846		Stand for body fluid-boxes	814034
			Small basket	814821
Adjustable handles	814001		Oxygen bottle holder	814820

Singapore (HQ)

Systmz Pte Ltd
No 1, Bukit Batok Crescent
WCEGA Plaza #09-07
Singapore 658064
Tel: +65 6465 0387
Fax: +65 6465 0656
Email: info@systmz.sg

Malaysia

Systmz Sdn Bhd
No.11, Jalan PJU 3/46
Sunway Damansara Technology Park
47810 Petaling Jaya
Selangor Darul Ehsan
Malaysia
Customer Care: 1300 222 500
Email: info@systmz.com.my

Dubai Liason Office

Al Leddan General Trading
P.O. Box 488176,
Dubai - UAE, Barsha
Dusseldorf Business Tower - 402
Tel: +971 4 4570375
Mobile: +971 50 2323001
Email: gm@leddan-ics.com